magnetic

M A L L E T

Medical device for surgery and dentistry

Instructions for use

TABLE OF CONTENTS

1.0. Important information for the user	3
2.0. Device Illustrations	4
2.1. Control elements, indicators, connections and functions	6
3.0. Explanation of the used symbols	7
4.0. Device description	9
4.1. Device general description	9
4.2. Description of the coupling parts	10
5.0. Safety Guidelines	12
5.1. Warning and Caution	12
5.2. Intended use	14
5.3. User's qualifications	15
5.4. Use conditions	15
5.5. Safety guidelines for the installation environment	15
5.6. Operating safety guidelines	15
6.0. Transport and storage	16
6.1. Transport and storage conditions	16
7.0. Installation and use	17
7.1. Unpacking	17
7.2. Installation and connection	18
7.3. Activation	19
7.4. Surgery instrument insertion on the handpiece OSTEOTOME function	19
7.5. Release of surgery instrument	20
7.6. How to activate the handpiece	21
7.7. Calibration of the force	21
8.0. Maintenance	22
8.1. Cleaning and maintenance	22
8.2. Disinfection and sterilization of the parts	22
8.3. Maintenance	24
8.4. Repairs	24
8.5. Fuses replacement	24
9.0. Disposal of waste	25
9.1. Responsibility	25
9.2. Warranty	25
9.3. Technical Literature	25
10.0. Technical description	26
10.1. Failures and malfunctions	26
11.0. Technical data	28
12.0. Electromagnetic compliance declaration	29

1.0. IMPORTANT INFORMATION FOR THE USER

This operating manual was made to help you to properly install and connect the Magnetic Mallet device. All the useful details for a proper use of this device are contained in this manual.

You should read it very carefully and store it in its slipcase in a dry and clean place in order to gather any useful information in the future.

Attention:

Please be aware you are referring to the last update of the instructions for use as the manufacturer reserves the right to make improvements at any time without prior warning. The images and drawings in this instructions for use are provided for illustrative purposes only.

© copyright Meta Ergonomica S.r.l.


All rights reserved. Meta Ergonomica S.r.l. retains all copyright in any text, graphic images and images owned by Meta Ergonomica S.r.l.

You may not mirror, modify or otherwise alter any file for re-broad, print or exploit the information contained therein without the written permission of Meta Ergonomica S.r.l.


2.0. DEVICE ILLUSTRATIONS

MAGNETIC MALLET ⁽¹⁾

(1) Magnetic Mallet can be also supplied with the Crown Remover Handpiece (on request)


Front view

Rear view

Osteotome Handpiece

2.1. Control elements, indicators, connections and functions

CONTROL UNIT

- A. Socket
- B. Mains plug
- C. ON/OFF switch
- D. Control unit
- E. Manual control push button
- F. Knob for the calibration of the force
- G. Display showing the force level
- H. Power cable
- I. Fuse holder
- L. Footswitch
- M. Handpiece connector
- N. Footswitch connector

OSTEOTOME HANDPIECE

- O. Osteotome function handpiece
- P. Chuck
- Q. Ring Nut

3.0. EXPLANATION OF THE USED SYMBOLS - LABELLING

SYMBOLS PLACED ON THE MAGNETIC MALLET CHASSIS

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
IN THE EVENT OF FAILURES CONTACT THE MANUFACTURER

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

CE 1936

SYMBOLS LABELLING

-  Information on the manufacturer and date of Company
-  Device serial number
-  AC (Alternating current)
-  Coupled part "BF type" (according to IEC 60601-1 rules)
-  Prior to setting the device going carefully read the instructions for use
-  Caution hazardous voltage
-  Store in a dry and clean place
-  Caution! The improper use of this device can cause injury
-  You should not use the device if the packaging is damaged
-  Storage temperature
-  CE mark including the identification number of the Notified Body
-  Disposal of special waste (electric and electronic devices)
-  Device Class II - IEC 417-5172

4.0. DEVICE DESCRIPTION

4.1. Device general description

Magnetic Mallet is an electro-medical device for surgery and dentistry practice with a double function according to the coupled part (handpiece O) :

- **OSTEOTOME** version to obtain the bone plastic deformation thanks to the available forces and their application timing.
- **CROWN REMOVER** version to obtain the crown and bridge removal from the abutment or the implant.

This device is basically made up by a control unit, a footswitch and two interchangeable handpieces which can be used depending on the needed practice.

The control unit case is made by self-extinguishing plastic. On the control unit you can find:

- The ON/OFF switch (C)
- Two different connectors. One for the footswitch and the other one for the handpiece. The two connectors are identifiable by the following graphic symbols:
 - ▶ Handpiece metal connector (M)
 - ▶ Footswitch plastic connector (N)
- Manual control push-button (E)
- Plastic knob to calibrate the level of the force (F)
- Display showing the selected force (G)

The control unit, the control/operational electronics together with the electronic/electric safeties are housed into the control unit case.

The handpiece is powered by the control unit with a 58V c.c. maximum voltage.

The handpiece is connected to the control unit by a cordset which can be manually unplugged. The handpiece can be activated by the pressure of both footswitch (L) and manual control push-button (E).


Both for the **OSTEOTOME** and for the **CROWN REMOVER** function the force applied to the handpiece through the pressure of the footswitch or of the push button is not continuous. Each pressure of the switch conveys a single force pulse

Turning the knob (F) you can select the desired level of the force among the four possible options (1-4); at the same time the selected force will be showed on the liquid crystal Display (G) and by the number selected on the knob.

Magnetic Mallet has been designed to give 1 pulses per second. You should not exceed the maximum threshold of 200 consecutive pulses without having a 40 minutes break.

4.2. Description of the coupling parts

Osteotome Function Handpiece with ring nut

For the OSTEOTOME function you must use the Hanpiece with the ring nut (O) equipped with the appropriate surgery instruments.

INTENDED USE

Medical device to perform a wide array of dentistry surgeries. The pulse conveyed to the surgery instrument allows you to achieve the plastic deformation of the bone in order to create the site for the implant.

The **Magnetic Mallet**, equipped with the handpiece and the surgery instrument is a medical device fitted to plastically deform both maxilla and lower jaw spongy-bone.

In this way the surgeon can create the proper site for the implant.

Through the handpiece the surgery instrument can be energized with 4 different forces so that the surgeon can replace the energy carried out by the human force in the manual practice (state of art).

You can select the desired level of the force among the 4 possible options (1-4) by turning the knob (F) placed on the control unit.

In the pre-marketing tests we carried out the most used forces have been the 2nd and the 3rd one; nevertheless it's important to remember that the whole spectrum of the forces has a very close progressive range. That's why they can all be considered suitable for surgery. It's entirely up to the surgeon to determine the more suitable one to achieve the desired result.

In order to familiarize yourself with the device, in the first stage you should always start from the lighter force (No.1). It's also advisable to keep the surgery instrument between your fingers selecting the different levels of the forces. In this way you shall be able to assess the corresponding pulse.

Surgery instrument action is led both by the pressure exerted by the surgeon on the handpiece and by the axial forwards of the energy. In any case the maximum forwards of the surgery instrument is of about 1 mm.

Attention: The surgeon must always determine the pressure to be exerted on the bone in order to achieve the desired result.

5.0. SAFETY GUIDELINES

5.1. Warning and Caution

Please carefully read this Manual for use and rigorously adhere to the instructions. The captions “Caution”, “Warning” and “Attention” have a precise meaning. You should carefully read the related instructions in order to grant a safe and effective working of this device.

Caution: indicates a danger/risk for the patient or the surgeon. The failure in comply with the indication can carry injuries to the patient or to the surgeon.

Warning: indicates the maintenance to comply with in order to avoid any possible damaging in the appliance.

Attention: indicates special information related to the use of the Magnetic Mallet device or other important information.

Caution: Prior to activating the device carefully read the manual for use. Pay close attention to the sections related to the device working and to safety guidelines in order to avoid any risk for the patient, your staff and yourself.

Caution: Install the Magnetic Mallet device keeping it out of the reach of the patient.

Caution: Always check the cables. If the cable is somehow damaged you should immediately replace it calling the authorized service provider.

Caution: Prior to activating the device you should always have a functional operational test (FOT) pressing the footswitch or the manual control push-button.

Caution: Prior to activating the device make sure you have connected the proper handpiece according to the intended practice/treatment.

Caution: It's strictly forbidden to open the control unit – Electric shock risk!
Any opening of the control unit by non authorized persons nullifies the warranty.

Caution: You should move away the device from the patient every time there is stalling in the device (for example: power failure during the use).

Caution: You should unplug the device from the main power every time you goes on working on the control unit (for example : disinfection/cleaning, unplugging of the connector).

Caution: It's forbidden to put fingers or any tool inside the footswitch or the hand-piece connectors.

Caution: Handle with care the packaging and avoid all sort of strikes or damages.

Caution: 1 pulse/sec for a maximum of 200 consecutive pulses then a 40 minutes break.

Warning: Keep dry! Make sure no liquid gets into the control unit. Do not lean liquid holders on it.

Warning: You should install the device in dry environments only and into the handpiece keep it dry. Take care no liquid gets into the control unit since this can lead to malfunction.

Warning: The device must be only supplied with the voltage which is stated on the rating plate placed on the rear panel of the control unit.

Warning: Use fuses having a value corresponding to the one stated on the line filter rating plate.

Warning: Handle and use with care the power cable and the hand-piece/footswitch cordsets. In order to avoid cable stress use them in wide leeway situations.

Attention: The warranty coverage is nullified by improper usage of the device.

Attention: According to an official test the Magnetic Mallet meets the requirements of Directive 2014/30/UE concerning electromagnetic compatibility and doesn't emit any interfering radiation.

Attention: Do not use on patients:

- Less than 14 years old
- With health problems
- With systemic chronic dependency
- With problems in maxillary sinus
- With handicaps
- With coagulation disorders
- With periodontal infection in alveolar bone
- Smokers
- Alcohol or drugs addicted

Attention: In the first stage you should start by a low level force and increase it step by step as necessary.

Attention: It's up to the surgeon to determine the more suitable pressure to be applied on the bone in order to achieve the desired result.

Attention: It's up to the surgeon to determine the more suitable force to achieve the desired result.

ATTENTION: You should be knowledgeable about the device and its working before using it on a patient.

5.2. Intended use

The **Magnetic Mallet** is a medical device to be used in dental surgery. According to the coupled handpiece it can have a different use:

- OSTEOTOME function HANDPIECE with ring nut
- CROWN REMOVER function* HANDPIECE with joint

Intended use : Osteotome function

Medical device to be used to perform dentistry surgeries. The pulse action plastically deforms the bone in order to implant partial denture.

The **Magnetic Mallet in OSTEOTOME function must be used with surgery instruments supplied by Meta Ergonomica only.**

* N.B. For the Crown-Remover function you should refer to the dedicated Instruction for use you can download from our website

For your own safety you should not use the device in other areas of applicability.

Attention: For your own safety you should not carry out alterations on the device.

5.3. User's qualifications

The Magnetic Mallet should be used in medical environments only and performed by a dentistry surgeon cognizant on the device use.

5.4. Use conditions:

Temperature -5°/+45°

5.5. Safety guidelines for the installation environment

The device should be used in dentist's surgery or medical environments only.

Warning: The device should be installed in dry environments only and kept dry.

Warning: The device must be only supplied with the voltage which is stated on the rating plate placed on the rear panel of the control unit.

5.6. Operating safety guidelines

Prior to your first treatment on a patient you should be cognizant on the working of this medical device.

Warning: 1 pulse/sec for a maximum of 200 consecutive pulses then a 40 minutes break

Warning: Prior to activating the device you should always carry a functional operational test (F.O.T.)

Warning: Instruments are not provided in a sterile package and must, therefore, be sterilised before their first use (see "Disinfection and sterilization of the parts" section). Carefully read the specific Sterilization handbook.

Warning: Move away the device from the patient every time there is a stalling in the device working (for example: power failure)

Warning: Prior to activate the device make sure you have plugged into the device the proper handpiece.

Warning: Prior to activate the device make sure you corerctly plugged the handpiece into the corresponding plug (M) and the footswitch in the corresponding plug (N).

Warning: Handle with care the blades supplied with the beaver. Blades by nature are very sharp-edged. You should pay attention in inserting them into the beaver. Do not insert beaver, with the blade coupled to it, into the handpiece, till the starting of the surgery. Do not leave the handpiece coupled with the beaver and the blade into the housing of the power supply.

Warning: make sure no liquids penetrate into the handpiece.

6.0. TRANSPORT AND STORAGE

6.1. Transport and storage conditions

Caution: Handle with care. Take care the package is not damaged.

Take care the external labelling is undamaged.

On the delivery please inspect the device into the packaging. Verify the device is undamaged and the wholeness of the supply (see “Unpacking” section).

7.0. INSTALLATION AND INSTRUCTIONS FOR USE

7.1. Unpacking

Take care in pulling the device and the accessories out from the package. Verify no article is missing or damaged.

Should any article be missing or be damaged please immediately contact the manufacturer or the supplier.

STANDARD SUPPLY

- 1 Control unit
- 1 Handpiece with ring nut – OSTEOTOME function
- 1 Footswitch
- 10 Surgery instruments
- 1 Spare nut
- 2 Spare O-rings for the ring nut
- 1 Instructions for use and maintenance

SUPPLY ONLY ON REQUEST

For further instruments download the catalogue from our website:
www.osseotouch.com

7.2. Installation and connection

Warning: Install the device in dry environments only. Take care no liquid gets into the control unit.

Caution: The device should be installed and kept away from the patient.

The power source and the Magnetic Mallet should have the same nominal voltage. The nominal voltage is indicated on the rating panel placed on the rear panel of the control unit.

2. Connect the power cable (H) to the device through the socket (A)
3. Plug in the power cable (H) of the device
4. Jack in the proper handpiece (O) into the connector (M)
5. Jack in the footswitch (L), into the connector (N)

Caution: Prior to the activation, take sure you've plugged in the proper handpiece according to the required use/treatment.


7.3. Activation

1. Switch on the device pushing the ON/OFF push button (C) placed on the rear panel of the control unit (D)
2. The display (G), placed on the front panel of the control unit, lights up.

Attention: When the display (G) lights up it will show a number (from 1 to 4) according to the positioning of the knob (F)

7.4 Surgery instrument insertion into the handpiece


1. Insert the instrument (R) into the handpiece
2. Insert the ring nut (Q) on the instrument
3. Screw the ring nut (Q)
4. Ensure the surgery instruments is tightly fastened in order to avoid any possible ejection of the instrument during surgery. This could cause severe risk for the operator, the assistant and the patient.

Warning: Handle with care the blades supplied with the beaver. Blades by nature are very sharp-edged. You should pay attention in inserting them into the beaver. Do not insert beaver, with the blade coupled to it, into the handpiece, till the starting of the surgery. Do not leave the handpiece coupled with the beaver and the blade into the housing of the power supply.

7.5. Release of surgery instrument

1. Unscrew the ring nut (Q)
2. Pull out the instrument (R)

7.6. How to activate the handpiece

1. To activate the handpiece, press the footswitch (L) or the manual control (E) which is placed on the front panel of the control unit (D).

Each pressure of the footswitch (L) or of the manual control (E) conveys a single pulse.

Attention: In the first stage you should always start by the lighter force (No.1).

Caution: 1 pulse/sec for a maximum of 200 consecutive pulses then a 40 minutes break

7.7. Calibration of the forces

1. Turn the knob (F) on the control unit front panel in order to calibrate the force level (from 1 to 4).

2. The display (G) shows the calibrated force level according the knob (F) positioning.

Attention: The most suitable level of the force is determined by the surgeon according to the practice/treatment.

The Surgery instrument action is mainly due to the pressure exerted by the surgeon on the handpiece and by the axial forwards movement of the energy.

Attention: It's up to the surgeon to determine the pressure to be exerted on the bone in order to achieve the desired result.

8.0. MAINTENANCE

8.1. Cleaning and maintenance

Caution: You should unplug the device from the main power every time you go working on the control unit (for example: cleaning, connectors unplugging).

Caution: The silicone cable with the connector and the handpiece should be sterilized by autoclave before any usage

Caution: You should never put your fingers or any tool inside the connector of the handpiece or into the footswitch one.

Warning: Take sure no liquids enter into the control unit. You should not hold liquids on it.

The control unit case should be routinely disinfected with cleaning products which aren't harmful for plastics.

Make sure no liquid enter into the handpiece

8.2. Disinfection and sterilization of the parts

Caution: You should unplug the device from the main power every time you go working on the control unit (for example: cleaning, connectors unplugging).

DISINFECTION AND STERILIZATION BY AUTOCLAVE OF THE HANDPIECE WITH ITS CORDSET AND CONNECTOR, THE RING NUT AND THE SURGERY INSTRUMENT.

The handpiece (O) with its cordset and metal connector (M), the ring nut (Q) and the surgery instrument (R) **must be submitted to autoclave sterilization before their first use and before all treatment.**

After having uncoupled the parts (see the section 7.5) you should carry out the cleaning process. Brush parts off under flowing tap water, then proceed with the sterilization by autoclave.

AUTOCLAVE DIRECTIONS FOR USE:

Achievement of the dew point through the proper combination of temperature and pressure values which should be the following ones:

- 121°C for 1 bar pressure
- 135°C for 2 bar pressure

ITEM	WASHING / DISINFECTION	STERILIZATION BY AUTOCLAVE
Handpiece with cordset and connector	Avoid anyt entering of liquids into the handpiece. Do not use liquid disinfectants for the handpiece	Items should be sterilized by autoclave before their first use as they're non-sterilized supplied. Items should be sterilized by autoclave before all treatment.
Chuck, joint	Prior to sterilization by autoclave you should disinfect the chuck and the joint. Soak a cotton bud with an hospital grade disinfectant and wipe any visible debris from both inside and outside the piece. To get a proper disinfection of the parts you should use hospital grade disinfectants. You should carefully read the instructions for use on the disinfectant packaging in order to determine the suitability of the treatment.	Items should be sterilized by autoclave before their first use as they're non-sterilized supplied. Items should be sterilized by autoclave before all treatment.
Osteotome/bone expanders, crown remover hook, inserter, cutter, beaver holder and beaver	Carry out the cleaning process. Brush off under flowing tap water then proceed with the sterilization by autoclave. Avoid any use of IUPAC (benzalkonium chloride) to clean instruments	Items should be sterilized by autoclave before their first use as they're non-sterilized supplied. Items should be sterilized by autoclave before all treatment.

Caution:

Prior to disinfection and sterilization, carefully read the instructions contained in the instruments tray.

8.3. Maintenance

Caution: Do not open the device! Electric shock danger! Any intervention on the device by unauthorized people nullifies the warranty.

NO LUBRICATION AND/OR MAINTENANCE IS NEEDED ON THIS DEVICE
NO INTERVENTION ON THIS DEVICE IS ALLOWED.

In the event of any failure you should immediately contact the manufacturer/distributor (see section “Failures and malfunctions”). The only maintenance you can carry out is the fuses replacement (see section “Fuses replacement”).

8.4. Repairs

Caution: Do not open the device! Electric shock risk! Any opening of this device by unqualified people nullifies the warranty.

Warning: FOR YOUR OWN SAFETY YOU SHOULD NOT MAKE ALTERATIONS TO THIS DEVICE.

IN THE EVENT OF FAILURES OR MALFUNCTIONS CALLING FOR REPAIRS YOU SHOULD PROMPTLY ADVISE THE MANUFACTURER/DISTRIBUTOR SO THAT THE DEVICE WILL BE CALLED IN.

The distributor/manufacturer will carry out the failures test and the repair works.


8.5. Fuses replacement

Warning: You should replace fuses with some new one having the amperage value indicated in the rating plate on the fuse holder.

1. Take off the fuse holder (I) from the line filter using a little screwdriver.
2. In the fuse holder you will find two 1,6 A fuses.
3. Verify if one or both fuses are burnt out.
4. Pull out the burnt out fuse/fuses.
5. Pull into the holder the new fuse/fuses. Make sure they have the same amperage value indicated in the rating plate placed on the fuse holder.
6. Verify the correct working of the device.

9.0. DISPOSAL OF WASTE

This device meets the requirements of **Directive 2012/19/UE** for electric and electronic devices disposal - **RAEE**. At the end of the period of use the device and the accessories, must be sent for recycling of the materials or for disposal in a manner which poses no threat to humans or the environment. The manufacturer is responsible for the compliance with National requirements.

9.1. Responsibility

The manufacturer is liable for safety, reliability and performances of this device only if:

- The installing has been performed in order to adhere rigorously to the instructions.
- All necessary alterations or repairs have been carried out by authorized repair service only.
- The device has been used in accordance with the instructions for use and its intended use.

9.2. Warranty

The manufacturer undertakes to provide the final customer of this device with a warranty of satisfactory functions freedom from faults in both and manufacturing process for the duration of 24 months from the delivery date. In case of justifiable complaints the manufacturer will provide repairs and/or spare parts free of charge. Nevertheless, the manufacturer will charge the final customer with shipment costs and it is not accountable for risks arising from the shipment itself. For other instances the manufacturer will refer to the warranty indicated in the trade general conditions.

Any opening, repair or alteration carried out by unauthorized persons relieves the manufacturer of all responsibility concerning the safe working of the device and nullifies the warranty.

9.3. Technical Literature

The manufacturer will furnish on request circuit diagram, the components list, all descriptions and information useful to the technical assistance in order to carry the authorized repairs.

10.0. TECHNICAL DESCRIPTION

Failures and malfunctions

Caution: Every time you go on working on the control unit you should unplug the device from the main power.

Failure / Malfunction

Break down without sonic alarm.

Possible cause

- Failure of electrical supply.
- Line filter fuse burnt out.
- Internal circuit failure.

Possible solution

- Verify the electrical supply presence.
- Fuse/fuses replacement (see section “Maintenance”). If the malfunction is still present after the fuse/fuses and power cable replacement contact the authorized service provider.
- Contact the authorized service provider.

Failure / Malfunction

The display shows the capital letter “E”, the sonic alarm is ringing and you notice a mild heating of the handpiece.

Possible cause

- Internal circuit failure.

Possible solution

- Switch off the device.
- Contact the authorized service provider.

Failure / Malfunction

When you press the footswitch no pulse is conveyed to the bone expander/osteotome.

Possible cause

- Possible wrong connection of the footswitch in the corresponding plug.
- The footswitch is damaged.
- Internal circuit failure.

Possible solution

- Verify you correctly connected the footswitch in the plug.
- Verify if the pulse is conveyed to the bone expander/osteotome by pressing the manual control push-button. If so contact the authorized service provider and ask for a new footswitch.
- Even if the footswitch doesn't work you can follow through the treatment using the manual control push-button.

Failure / Malfunction

If you press the manual control push-button no pulse is conveyed to the bone-expander/osteotome.

Possible cause

- The manual control push-button is damaged.
- Internal circuit failure

Possible solution

- Verify if the pulse is conveyed by pressing the footswitch. If so contact the authorized service asking for the manual control switch repair/replacement.
- Contact the authorized service.

Failure / Malfunction

The display is switched off or it shows some nondescript number.

Possible cause

- The display is damaged.
- Internal circuit failure.

Possible solution

- Take note that the level of the force is all the same indicated by the knob positioning (1-2-3-4). In that event you can follow through the treatment before calling the authorized service for display replacement or repairs.

Failure / Malfunction

Mild heating of the handpiece

Possible cause

- Improper use. You have exceeded the maximum threshold of 200 consecutive pulses without having had the necessary break.

Possible solution

- Switch off the device. Wait 40 minutes before switching on the device again.

Failure / Malfunction

The o-ring type 108 of the handpiece is no more working as anti-unscrewing for the ring nut

Possible cause

- O-ring consumption or breaking

Possible solution:

- Replace the O-ring with the spare one supplied with the handpiece (see blades container).
- If the problem endures, contact the authorized service provider.

11.0. TECHNICAL DESCRIPTION

Supply voltage:	230 V
Frequency:	50 Hz
Current input (stand by):	0,070 A
Current input (operative condition):	0,38 A
Fuse:	F - 5x20 - 250V - 1.6A

Medical device class IIa in accordance with MDD 93/42/CEE

Medical device Class II – BF type in accordance with IEC 60601-1

Dimensions (l x h x d)	226.5 mm x 223 mm x 109 mm
Control unit weight	Kg. 2,600
Handpiece weight	Kg. 0,390

Transport/storage conditions

Storage temperature range	- 40°C / + 50°
Humidity (relative humidity, without condensation)	10% . 90%
Air pressure:	+700hPa / +1060 hPa

Conditions for use

Temperature range for use	-5°C / + 45°C
Humidity (relative humidity, without condensation)	15% . 85%
Air pressure:	+700hPa / +1060 hPa

12.0. ELECTROMAGNETIC COMPLIANCE DECLARATION

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, 4, 8, 15kV air	Compliant	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/Burst IEC 61000-4-4	±2kV for input power ports ±1kV for I/O ports	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Input power ports: 0.5 and 1.0kV (line to line) 0.5, 1.0 and 2.0kV (line to earth) Signal I/O: 2kV (line to earth)	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Dips: >95% UT for 0.5 and 1 cycle 30% UT for 25 cycles (50Hz) 30% UT for 30 cycles (60 Hz) Interruption: >95% UT for 250 cycles (50Hz) >95% UT for 300 cycles (60Hz)	Compliant	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions it is recommended that the device be powered from and uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Compliant	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150KHz to 80MHz</p> <p>3V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> <p>$d=1,167* \text{sqrt} (P)$</p> <p>$d=1,167* \text{sqrt} (P)$ 80 MHz to 800 MHz $d=2,2333* \text{sqrt} (P)$ 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey³, should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: at 80 MHz and 800 MHz, the higher frequency range applies

Note 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The equipment is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum power of communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 KHz to 80 MHz $d=1,17*\sqrt{P}$ m	80 MHz to 800 MHz $d=1,17*\sqrt{P}$ m	800 MHz to 2,5 GHz $d=2,33*\sqrt{P}$ m
0,01	0,117	0,117	0,233
0,1	0,370	0,370	0,740
1	1,17	1,17	2,33
10	3,70	3,70	7,40
100	11,7	11,7	23,3

For transmitters rated at maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150KHz to 80MHz</p> <p>3V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> <p>$d=1,167* \text{sqrt}(P)$</p> <p>$d=1,167* \text{sqrt}(P)$ 80 MHz to 800 MHz $d=2,2333* \text{sqrt}(P)$ 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey³, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

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Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, 4, 8, 15kV air	Compliant	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/Burst IEC 61000-4-4	±2kV for input power ports ±1kV for I/O ports	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Input power ports: 0.5 and 1.0kV (line to line) 0.5, 1.0 and 2.0kV (line to earth) Signal I/O: 2kV (line to earth)	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Dips: >95% UT for 0.5 and 1 cycle 30% UT for 25 cycles (50Hz) 30% UT for 30 cycles (60 Hz) Interruption: >95% UT for 250 cycles (50Hz) >95% UT for 300 cycles (60Hz)	Compliant	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Compliant	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150KHz to 80MHz</p> <p>3V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> <p>$d=1,167* \text{sqrt} (P)$</p> <p>$d=1,167* \text{sqrt} (P)$ 80 MHz to 800 MHz $d=2,2333* \text{sqrt} (P)$ 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey³, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

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© Copyright Meta Ergonomica S.r.l.
MANUFACTURER: META ERGONOMICA S.r.l.